



AMCHOMATRIX

STERILE DEHYDRATED HUMAN
AMNION - INTERMEDIATE LAYER -
CHORION MEMBRANE ALLOGRAFT



OVERVIEW

AmchoMatrix can be used for homologous use as a wound covering for non-infected, acute and chronic wounds that occur due to certain conditions

CELLUTION BIOLOGICS

 (888)575-7357

 www.cellutionbiologics.com

TREATS WOUNDS AS RESULT OF:

- Diabetes
- Burns
- Peripheral Vascular Arterial Disease
- Post-Operative Wounds
- Chronic Venous Insufficiency Post Traumatic Wounds

ABOUT

AmchoMatrix is created using placental membranes that have significant applications in fields such as wound healing, tissue repair, and regenerative medicine. The safety and quality of these biological materials is ensured by obtaining informed consent, thorough health screening, and physical examinations to minimize risks for both donors and AmchoMatrix recipients.

USAGE GUIDELINES

- Wound site should be properly prepared for debridement.
- The wound area should be cleared from any tissue or infection that may hinder wound healing.
- The AmchoMatrix wound graft is applied with notch facing towards the wound bed and secured with non-adhesive dressings.
- Available in 2x3 cm, 3x3 cm, 4x4 cm, 4x6 cm, 4x8 cm, and 10x10 cm.